


## T MAZE



Note: Drawing used for illustration purposes.

### Features

- Made of durable and easy to clean ABS plastic
- Available in black, white or beige
- One surface of the maze is smooth; the other is slightly textured
- Designed for optimal use with SDI video tracking systems that feature:
  - User defined zones
  - Tracks time spent and distance traveled, speed and resting time in each zone
  - Customizable Reports

### T MAZE SPECIFICATIONS

|                      | BASE                      | MAZE DIMENSIONS |        |         |          |        |
|----------------------|---------------------------|-----------------|--------|---------|----------|--------|
| Mouse                | 22" x 17.5"               | (A) 14"         | (B) 2" | (C) 6"  | (D) 7.5" | (E) 5" |
| Rat                  | 36" x 27"                 | (A) 28"         | (B) 4" | (C) 12" | (D) 15"  | (E) 8" |
| Material Composition | 1/4" Haircell ABS Plastic |                 |        |         |          |        |
| Color Options        | Black, White, Beige       |                 |        |         |          |        |

| PRODUCT NAME    | PART NUMBER     | PRODUCT NAME    | PART NUMBER     |
|-----------------|-----------------|-----------------|-----------------|
| T Maze for Mice | Black 7001-0236 | T Maze for Rats | Black 7001-0237 |
|                 | White 7001-0329 |                 | White 7001-0330 |
|                 | Beige 7001-0331 |                 | Beige 7001-0332 |

NOTE: Beige is SDI's standard Maze color unless specified by the researcher. The SMART Video Tracking System Software is sensitive enough to track a subject with little or no contrast.