


Y MAZE



Features & Benefits

- Made of durable and easy to clean ABS plastic
- Available in black, white or beige
- One surface of the maze is smooth; the other is slightly textured
- Designed for optimal use with SDI video tracking systems that feature:
 - User defined zones
 - Tracks time spent and distance traveled, speed and resting time in each zone
 - Customizable Reports

Note: Drawing used for illustration purposes.

Y MAZE SPECIFICATIONS

	BASE	MAZE DIMENSIONS			
Mouse	18" x 18"	(A) 2"	(B) 6"	(C) 8"	(D) 5"
Rat	30" x 30"	(A) 4"	(B) 12"	(C) 16"	(D) 8"
Material Composition	1/4" Haircell ABS Plastic				
Color Options	Black, White, Beige				

PRODUCT NAME	PART NUMBER	PRODUCT NAME	PART NUMBER
Y Maze for Mice	Black 7001-0307	Y Maze for Rats	Black 7001-0304
	White 7001-0308		White 7001-0305
	Beige 7001-0309		Beige 7001-0306

NOTE: Beige is SDI's standard Maze color unless specified by the researcher. The SMART Video Tracking System Software is sensitive enough to track a subject with little or no contrast.