



OPEN FIELD SINGLE UNIT

PART NUMBER

BEIGE


➤ 7001-0357

BLACK

➤ 7001-0355

WHITE

➤ 7001-0356



Drawings and images are not shown proportionally for illustration purposes.

FEATURES & BENEFITS

Made of durable and easy to clean ABS plastic

The interior walls of the maze are smooth, while the floor has a slightly textured surface

Walls and floors separate for easy cleaning

Designed for optimal use with SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports

Open Field Single Unit Specifications

Dimensions	
Inside Wall Length (A)	20"
Inside Chamber Footprint	20" x 20"
Inside Wall Height (B)	15.25"
Base of Maze	22" x 22"
Material Composition	.25" ABS Plastic
Color Options	Black, White, Beige

For more information on any of our products or services please visit us on the Web at:

www.sandiegoinstruments.com

or contact us via email at

sales@sandiegoinstruments.com



9155 Brown Deer Rd Suite 8
San Diego CA 92121
Phone 858-530-2600
Fax 858-530-2626