



# BARNES MAZE

## MICE PART NUMBER

### BEIGE

➤ 7001-0340

### BLACK

➤ 7001-0324

### WHITE

➤ 7001-0235

## RAT PART NUMBER

### BEIGE

➤ 7001-0341

### BLACK

➤ 7001-0325

### WHITE

➤ 7001-0326



## FEATURES & BENEFITS

Made of durable and easy to clean ABS plastic

One target box and 20 false boxes

Available in black, white or beige

Designed for optimal use with SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports

## Barnes Maze Specifications

Dimensions	
Mouse Diameter	36"
Rat Diameter	48"
Mouse No. of Holes	Qty 20 2" Diam.
Rat No. of Holes	Qty 20 4" Diam.
Mouse Start Tube	8" (H) 6" (D)
Rat Start Tube	12" (H) 8" (D)
Material Composition	
Top	.25" ABS
Stand	Cast Iron
Boxes	.125" PVC Foam
Color Options	White, Beige or Black
Weight	
Rat	48lbs
Mouse	25 lbs

For more information on any of our products or services please visit us on the Web at:

[www.sandiegoinstruments.com](http://www.sandiegoinstruments.com)

or contact us via email at

[sales@sandiegoinstruments.com](mailto:sales@sandiegoinstruments.com)

