



ELEVATED PLUS MAZE

MICE PART NUMBER

BEIGE

➤ 7001-0316

BLACK

➤ 7001-0225

WHITE

➤ 7001-0336

RAT PART NUMBER

BEIGE


➤ 7001-0337

BLACK

➤ 7001-0227

WHITE

➤ 7001-0338



Drawings and images are not shown proportionally for illustration purposes.

FEATURES & BENEFITS

Made of durable and easy to clean ABS plastic

One surface of the maze is smooth; the other is slightly textured

Available in black, white or beige

Designed for optimal use with SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports

Elevated Plus Maze Specifications

Dimensions	
Mouse	(A) 26" x (B) 12" x (C) 2" x (D) 15.25" (E) 21.25" (F) .125" (lip)
Rat	(A) 43.5" x (B) 19.75" x (C) 4" x (D) 19.5" x (E) 31.5" x (F) .35" (lip)
Material Composition	.25" ABS plastic
Color Options	White, Beige or Black

For more information on any of our products or services please visit us on the Web at:
www.sandiegoinstruments.com
or contact us via email at
sales@sandiegoinstruments.com



9155 Brown Deer Rd Suite 8
San Diego CA 92121
Phone 858-530-2600
Fax 858-530-2626