



# OPEN FIELD CONVERTIBLE

## FEATURES & BENEFITS


Made of durable and easy to clean ABS plastic

One surface of the maze is smooth; the other is slightly textured

Available in black, white or beige

Designed for optimal use with SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports



*Drawings and images are not shown proportionally for illustration purposes.*

## Open Field Convertible Specifications

Dimensions	
Base	43" x 42"
Outer Wall Length	40" x 40"
Wall Height	15"
Inside Chamber	19.5" x 19.5"
Material Composition	.25" ABS Plastic
Color Options	Black, White, Beige

### PART NUMBER

BEIGE

➤ 7001-0067

BLACK

➤ 7001-0353

WHITE

➤ 7001-0354

For more information on any of our products or services please visit us on the Web at:

[www.sandiegoinstruments.com](http://www.sandiegoinstruments.com)

or contact us via email at

[sales@sandiegoinstruments.com](mailto:sales@sandiegoinstruments.com)



9155 Brown Deer Rd Suite 8  
San Diego CA 92121  
Phone 858-530-2600  
Fax 858-530-2626