



RADIAL ARM MAZE

MICE PART NUMBER

BEIGE

➤ 7001-0443

BLACK

➤ 7001-0440

WHITE

➤ 7001-0437

RAT PART NUMBER

BEIGE


➤ 7001-0452

BLACK

➤ 7001-0449

WHITE

➤ 7001-0446



Drawings and images are not shown proportionally for illustration purposes.

FEATURES & BENEFITS

Made of durable and easy to clean ABS plastic

Interchangeable plain and pellet feeder floors

Includes manual doors at each arm

Available in black, white or beige

Designed for optimal use with SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports

Radial Arm Maze Specifications

Dimensions	
Mouse	(A) 26.45" x (B) 9" x (C) 8.45" x (D) 3" x (E) 6"
Rat	(A) 47.86" x (B) 18.5" (C) 10.86" x (D) 4" x (E) 10"
Mouse Pellet Feeder	(H) .75" Diam x .5" Deep
Rat Pellet Feeder	(H) 1" Diam x .75" Deep
Material Composition	.25" ABS Plastic
Color Options	Black, White, Beige

For more information on any of our products or services please visit us on the Web at: www.sandiegoinstruments.com or contact us via email at sales@sandiegoinstruments.com



9155 Brown Deer Rd Suite 8
San Diego CA 92121
Phone 858-530-2600
Fax 858-530-2626