



# SOCIABILITY APPARATUS


## PART NUMBER

SOCIABILITY APPARATUS

➤ 7500-0362

STRANGER CAGE

➤ 7500-0397



## FEATURES & BENEFITS

**Three chamber enclosure to allow for easy exploration and social behavior**

**Made of Acrylic walls with ABS plastic floor**

**Weighted wire cages provide isolation for novel animal**

**Removable walls allow for versatility using just one enclosure**

**Capture all movements with the ANY-maze video tracking system**

## Sociability Apparatus Specifications

<b>Dimensions</b>	
3 Compartments	24" (W) x 16" (D) x 9" (H)
Door Dimension	2" (W) x 3.5" (H)
Stranger Enclosure	2.75" (D) x 6.75" (H)
Stranger Enclosure Bars	.125" (D) x .25" (Spacing)
<b>Material Composition</b>	
Walls	Acrylic
Floor	ABS Plastic

For more information on any of our products or services please visit us on the Web at:

[www.sandiegoinstruments.com](http://www.sandiegoinstruments.com)

or contact us via email at

[sales@sandiegoinstruments.com](mailto:sales@sandiegoinstruments.com)



9155 Brown Deer Rd Suite 8  
San Diego CA 92121  
Phone 858-530-2600  
Fax 858-530-2626