



T MAZE

MICE PART NUMBER

BEIGE

➤ 7001-0407

BLACK

➤ 7001-0404

WHITE

➤ 7001-0401

RAT PART NUMBER

BEIGE

➤ 7001-0416

BLACK

➤ 7001-0413

WHITE

➤ 7001-0410

For more information on any of our products or services please visit us on the Web at:

www.sandiegoinstruments.com

or contact us via email at

sales@sandiegoinstruments.com

FEATURES & BENEFITS

Made of durable and easy to clean ABS plastic


Includes manual doors at each upper arm and starter door in base leg

Interchangeable plain and pellet feeder floors

Available in black, white or beige

Designed for optimal use with SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports



Drawings and images are not shown proportionally for illustration purposes.

T Maze Specifications

| Dimensions | |
|----------------------|--|
| Mouse | (A) 24" x (B) 3" x (C) 10.5" x (D) 15" x (E) 5" x (F) 3.5" |
| Rat | (A) 36" x (B) 4" x (C) 16" x (D) 22" x (E) 10" x (F) 7.25" |
| Mouse Pellet Feeder | (H) .75" Diam x .5" Deep |
| Rat Pellet Feeder | (H) 1" Diam x .75" Deep |
| Material Composition | .25" ABS Plastic |
| Color Options | Black, White, Beige |



9155 Brown Deer Rd Suite 8
San Diego CA 92121
Phone 858-530-2600
Fax 858-530-2626