



Y MAZE

MICE PART NUMBER

BEIGE

➤ 7001-0425

BLACK

➤ 7001-0422

WHITE

➤ 7001-0419

RAT PART NUMBER

BEIGE


➤ 7001-0434

BLACK

➤ 7001-0431

WHITE

➤ 7001-0419



Drawings and images are not shown proportionally for illustration purposes.

FEATURES & BENEFITS

Made of durable and easy to clean ABS plastic

Includes manual doors at each upper arm and starter door in base leg

Interchangeable plain and pellet feeder floors

Available in black, white or beige

Designed for optimal use with SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports

Y Maze Specifications

Dimensions	
Mouse	(A) 3" x (B) 15" x (C) 15" x (D) 5" x (E) 3.5"
Rat	(A) 4" x (B) 22" x (C) 22" x (D) 10" x (E) 7"
Mouse Pellet Feeder	(H) .75" Diam x .5" Deep
Rat Pellet Feeder	(H) 1" Diam x .75" Deep
Material Composition	.25" ABS Plastic
Color Options	Black, White, Beige

For more information on any of our products or services please visit us on the Web at:

www.sandiegoinstruments.com

or contact us via email at

sales@sandiegoinstruments.com



9155 Brown Deer Rd Suite 8
San Diego CA 92121
Phone 858-530-2600
Fax 858-530-2626