



# ZERO MAZE

## MICE PART NUMBER

### BEIGE

➤ 2325-0231

### BLACK

➤ 2325-0230

### WHITE

➤ 2325-0229

## RAT PART NUMBER

### BEIGE

➤ 2325-0228

### BLACK

➤ 2325-0227

### WHITE


➤ 2325-0226

For more information on any of our products or services please visit us on the Web at:

[www.sandiegoinstruments.com](http://www.sandiegoinstruments.com)

or contact us via email at

[sales@sandiegoinstruments.com](mailto:sales@sandiegoinstruments.com)



\*Drawings and images are not shown proportionally for illustration

## FEATURES & BENEFITS

Made of durable and easy to clean ABS plastic

One surface of the maze is smooth; the other is slightly textured

Available in black, white or beige

Designed for optimal use with SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports

## Zero Maze Specifications

Dimensions	
Mouse	(A) 20" x (B) 24" x (C) 6"
Rat	(A) 40" x (B) 48" x (C) 8"
Adjustable Legs	
Mouse	3 Legs 19.75" - 20.75"
Rat	6 Legs 19.75" - 20.75"
Material Composition	.25" ABS Plastic
Color Options	Black, White, Beige



9155 Brown Deer Rd Suite 8  
San Diego CA 92121  
Phone 858-530-2600  
Fax 858-530-2626