



BARNES MAZE

FEATURES & BENEFITS

Made of durable and easy to clean ABS plastic

Used for spatial learning and memory

One target box and 20 false boxes

Available in black, white or beige

SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports



MICE PART NUMBER

BEIGE

➤ 7001-0340

BLACK

➤ 7001-0324

WHITE

➤ 7001-0235

RAT PART NUMBER

BEIGE

➤ 7001-0341

BLACK

➤ 7001-0325

WHITE

➤ 7001-0326

Barnes Maze Specifications

Dimensions	
Mouse Diameter	36"
Rat Diameter	48"
Mouse No. of Holes	Qty 20 2" Diam.
Rat No. of Holes	Qty 20 4" Diam.
Mouse Start Tube	8" (H) 6" (D)
Rat Start Tube	12" (H) 8" (D)
Material Composition	
Top	.25" ABS
Stand	Cast Iron
Boxes	.125" PVC Foam
Color Options	White, Beige or Black
Weight	
Rat	48lbs
Mouse	25 lbs

For more information on any of our products or services please visit us on the Web at:

www.sandiegoinstruments.com

or contact us via email at

sales@sandiegoinstruments.com



9155 Brown Deer Rd Suite 8
San Diego CA 92121
Phone 858-530-2600
Fax 858-530-2626