



WATER TIGHT T MAZE

MICE PART NUMBER

BEIGE

➤ 7001-0484

BLACK

➤ 7001-0482


WHITE

➤ 7001-0483

ESCAPE PLATFORM

ACRYLIC

➤ 7001-0398



Drawings and images are not shown proportionally for illustration purposes.

FEATURES & BENEFITS

Made of durable and easy to clean ABS plastic

Includes manual doors at each upper arm and starter door in base leg

Glued plain floor and walls assure maze is watertight

Available in black, white or beige

Designed for optimal use with SDI ANY-maze video tracking system that features:

- User defined zones
- Tracks time spent and distance traveled, speed and resting time in each zone.
- Customizable Reports

Water Tight T Maze Specifications

Dimensions	
Mouse	(A) 24" x (B) 3" x (C) 10.5" x (D) 15" x (E) 12" x (F) 3.5"
Material Composition	.25" ABS Plastic
Color Options	Black, White, Beige
Accessories	
Escape Platform	3" (W) X 3" (D) 6" X (H)
Material Composition	Acrylic

For more information on any of our products or services please visit us on the Web at: www.sandiegoinstruments.com or contact us via email at sales@sandiegoinstruments.com



9155 Brown Deer Rd Suite 8
San Diego CA 92121
Phone 858-530-2600
Fax 858-530-2626